


www.falconmedical.eu

D₃

Dental Equipment



Falcon[®] [™]
Medical Polska
Your Instruments since 1975

Falcon *K-Series* 1 Year Guarantee

High speed air turbine fibre optic handpieces

- Ceramic bearings
- Quiet operation
- Perfect balance
- Anti-retraction valve
- Push button chuck
- Quick coupling fitting



DQ.2100.00



DQ.2100.25



Specifications

Code	Head	Head size	Maximum speed	Cooling	Bur size
DQ.2100.00	Torque	14.6 mm	300,000~ 380,000 R.P.M.	4-Spray	Ø 1.59 - 1.60 mm , Length Standard FG
DQ.2100.25	Mini	12.6 mm	300,000~ 380,000 R.P.M.	Single	Ø 1.59 - 1.60 mm , Length Standard FG

Quick swivel coupling



DQ.2525.00

Fibre optic system



DQ.2610.00

2 handpieces output unit Transformer

Silicone tubing for fibre optic handpieces



DQ.2600.06

2,8 meters with adapter



Falcon

N-Series

1 Year Guarantee

High speed air turbine handpieces



- Ceramic bearings
- Quiet operation
- Perfect balance



Fully autoclavable

Spare parts are available, please ask for details

Falcon N-Series

Torque head handpieces

Torque head

- Anti-retraction valve
- Ceramic bearings
- Push button chuck
- Quick coupling fitting



DQ.2200.00 Mini head



DQ.2225.00 Standard head

Specifications

Code	Head size	Length	Recommended air pressure	Maximum speed	Bur size	Length	Single-spray coolant
DQ.2200.00	∅ 11.2 mm	13.2mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000-420,000 R.P.M.	∅ 1.59 - 1.60 mm	Standard	✓
DQ.2225.00	∅ 13.2 mm	14.6mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	300,000-350,000 R.P.M.	∅ 1.59 - 1.60 mm	Standard	✓

Quick swivel couplings

- Anti-retraction valve
- 360° swivels
- Fully autoclavable



DQ.2500.02
2-holes Borden



DQ.2500.04
4-holes Midwest

Standard head handpieces

Standard Torque, push button chuck



DQ.2325.00 Standard head

- Anti-retraction valve
- Ceramic bearings
- Quick coupling fitting



Standard Torque, wrench type chuck



DQ.2375.00 Standard head



Standard Torque, push button chuck, ceramic bearing



DQ.2425.02 2-holes Borden
DQ.2425.04 4-holes Midwest
Standard Head

Standard Torque, wrench type chuck, ceramic bearing



DQ.2475.02 2-holes Borden
DQ.2475.04 4-holes Midwest
Standard Head

Specifications

Code	Head size	Length	Recommended air pressure	Maximum speed	Bur size	Length	Single- spray coolant
DQ.2325.00	∅ 11.2 mm	14.6mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000~ 420,000 R.P.M.	∅ 1.59 ~ 1.60 mm	Standard	✓
DQ.2375.00	∅ 10.2 mm	13.6mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000~ 420,000 R.P.M.	∅ 1.59 ~ 1.60 mm	Standard	✓
DQ.2425.02	∅ 11.2 mm	14.6mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000~ 420,000 R.P.M.	∅ 1.59 ~ 1.60 mm	Standard	✓
DQ.2425.04	∅ 11.2 mm	14.6mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000~ 420,000 R.P.M.	∅ 1.59 ~ 1.60 mm	Standard	✓
DQ.2475.02	∅ 10.2 mm	13.2 mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000~ 420,000 R.P.M.	∅ 1.59 ~ 1.60 mm	Standard	✓
DQ.2475.04	∅ 10.2 mm	13.2 mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	320,000~ 420,000 R.P.M.	∅ 1.59 ~ 1.60 mm	Standard	✓

Silicone handpiece tubing, 2-holes Borden



DQ.2600.02
2.8 meters with adapter

Silicone handpiece tubing, 4-holes Midwest



DQ.2600.04
2.8 meters with adapter

Contra angle handpiece 1:1, push button chuck

1 Year Guarantee



DQ.2800.00

Specifications

Code	Head size	Length	Recommended air pressure	Maximum speed	Bur size	Single- spray coolant
DQ.2800.00	∅ 8 mm	13.7 mm	2.0-2.5 kg / cm ² ; (0.2- 0.25 Mpa)	350,000 R.P.M.	∅ 2.35 mm	✓

3-1 Air water syringes and accessories

1 Year Guarantee



DQ.2705.00
Complete with tip



DQ.2700.00
Complete with 2 tips

3-1 Air water syringe tips



DQ.2707.00



DQ.2708.00
Disposable
Pack of 150 pieces

Polyurethane tube for air water syringe, 2-holes



DQ.2710.02
2 meters long



Air polisher handpieces

- For intraoral use
- Light weight and ergonomic design
- Detachable nozzle



Manual cleaning and cold sterilization only

1 Year Guarantee

P-Mate-I




DQ.3000.02 2-holes Borden
 DQ.3000.04 4-holes Midwest

Polishing powder for P-Mate handpieces



DQ.3001.01 Lemon
 DQ.3001.02 Mint
 DQ.3001.03 Strawberry
 130 g

P-Mate-II



DQ.3000.06 2-holes Borden
 DQ.3000.08 4-holes Midwest



Sealing machine and accessories

1 year guarantee for equipment, heating element or other parts may need replacement depending on frequency of use or misuse.

- Sealing width 10 mm
- Maximum width of sealing 200 mm



DQ.9000.01

Code	Product
DQ.9000.02	Roll stand only
DQ.9000.03	Worktable only
DQ.9000.04	Heating element
DQ.9000.05	Cutting blade and handle
DQ.9000.06	Teflon sealing tape
DQ.9000.07	Upper rubber strip

Microblaster for extraoral use

1 Year Guarantee



DQ.3005.01



DQ.9000.00
 Complete with roll stand and working table