


*D*₁₂

Impression Accessories





Impression Accessories

Unipress®



all for 1 – 1 for all

DI.1150.00

Universal dispenser kit for all ratios of 25ml and 50ml cartridges

Contents

- 1 1 x Dispenser
- 2 1 x Cartridge holder for rectangular flange cartridges (old type)
- 3 1 x Cartridge holder for round flange cartridges (new type)
- 4 1 x Plunger for 1:1/2:1 ratio
- 5 1 x Plunger for 4:1 ratio
- 6 1 x Plunger for 10:1 ratio



DENTALINE®

The Single-Use Product Specialists™



Mixing tips 1:1 ratio



- DI.1106.01 Green, pack of 50
- DI.1106.02 Pink, pack of 50
- DI.1106.03 Yellow, pack of 50
- DR.9815.01 Brown, pack of 50
- DR.9815.02 Brown, pack of 50
- DR.9814.01 Brown, pack of 50

For 50 ml round flange dual cartridges

For 25 ml round flange dual cartridges

Disposable impression syringes



DI.016.000
Pack of 25



- DI.1100.00 Clear, pack of 50
- DI.1100.01 Clear, pack of 50
- DI.1100.03 Green, pack of 50
- DI.1100.04 Pink, pack of 50
- DI.1100.02 Yellow, pack of 50
- DR.9801.00 Clear, pack of 50
- DR.9800.00 Clear, pack of 50


Bite gauge

A simple measuring instrument as an aid in the re-establishment of correct facial contour during full or partial denture work.



DI.001.000
Falcon

Slide **B** is fully adjustable and reversible along shaft **A** and may be fixed in any position by means of its locking-screw as shown.



Bite plane



DI.002.000
Fox

Carvers for silicone impression materials



DI.010.017




DI.010.014

Gingival cord packers




DI.005.010 serrated



DI.009.010 serrated



DI.005.020 smooth



DI.006.020 serrated

Gingival cord packers



DI.005.014 serrated



DI.005.024 smooth



Disposable impression trays

DENTALINE

Disposable impression trays, available in green or blue colors



Pack includes 25 trays of 1 size and 1 handle

Trial Pack

- DI.025.000
- DI.026.000

25 assorted trays and 1 handle

Edentulous



- | | | | |
|--|------------------------------------|------------------------------------|------------------------------------|
| Fig. 3
● DI.025.030
● DI.026.030 | Fig. 4
DI.025.040
DI.026.040 | Fig. 5
DI.025.050
DI.026.050 | Fig. 6
DI.025.060
DI.026.060 |
| Large | | Medium | |

Dentate



- | | | | |
|---|-------------------------------------|-------------------------------------|-------------------------------------|
| Fig. 11
● DI.025.110
● DI.026.110 | Fig. 12
DI.025.120
DI.026.120 | Fig. 13
DI.025.130
DI.026.130 | Fig. 14
DI.025.140
DI.026.140 |
| Large | | Medium | |

Orthodontic



- | | | | |
|---|-------------------------------------|-------------------------------------|-------------------------------------|
| Fig. 15
● DI.025.150
● DI.026.150 | Fig. 16
DI.025.160
DI.026.160 | Fig. 17
DI.025.170
DI.026.170 | Fig. 18
DI.025.180
DI.026.180 |
| Large | | Medium | |

Partial



- | | | |
|---|---|---|
| Fig. 20
● DI.025.200
● DI.026.200
Anterior | Fig. 21
DI.025.210
DI.026.210
Upper Left / Lower Right | Fig. 22
DI.025.220
DI.026.220
Upper Right / Lower Left |
|---|---|---|

Spare metal handle



DI.020.010
Pack of 10 pieces



DENTALINE

Disposable transparent impression trays
Pack of 10 pieces



- | | | | |
|---------------------|---------------------|---------------------|---------------------|
| Fig.1
DI.029.010 | Fig.2
DI.029.020 | Fig.3
DI.029.030 | Fig.4
DI.029.040 |
|---------------------|---------------------|---------------------|---------------------|

Trial Pack
DI.029.000
9 assorted trays



- | | | | | |
|---------------------|---------------------|---------------------|---------------------|---------------------|
| Fig.5
DI.029.050 | Fig.6
DI.029.060 | Fig.7
DI.029.070 | Fig.8
DI.029.080 | Fig.9
DI.029.090 |
|---------------------|---------------------|---------------------|---------------------|---------------------|

DENTALINE

Disposable impression trays

For simultaneous crown and bridge impression,
counter impression and bite registration.
Pack of 10 pieces

Trial Pack
DI.039.000
10 assorted trays



- | | | |
|---------------------|---------------------|---------------------|
| Fig.1
DI.039.010 | Fig.2
DI.039.020 | Fig.3
DI.039.030 |
|---------------------|---------------------|---------------------|



- | | |
|---------------------|---------------------|
| Fig.4
DI.039.040 | Fig.5
DI.039.050 |
|---------------------|---------------------|

1-2-3 impression trays

For simultaneous crown and bridge impression,
counter impression and bite registration.



- | | | | | |
|---------------------|---------------------|---------------------|---------------------|---------------------|
| Fig.1
DI.040.010 | Fig.2
DI.040.020 | Fig.3
DI.040.030 | Fig.4
DI.040.040 | Fig.5
DI.040.050 |
|---------------------|---------------------|---------------------|---------------------|---------------------|

D-Color Impression Trays

- Color Coded for easy identification of size
- Micro etched inner surface for better retention
- Available perforated or solid









ID-Color Regular Impression trays with rim, perforated

Lower

					
L1 / XL DI.801.010 Purple	L2 / L DI.801.020 Green	L3 / M DI.801.030 Blue	L4 / S DI.801.040 Red	L5 / XS DI.801.050 Yellow	DI.851.060 Set of 6 pieces

Upper

					
U1 / XL DI.811.010 Purple	U2 / L DI.811.020 Green	U3 / M DI.811.030 Blue	U4 / S DI.811.040 Red	U5 / XS DI.811.050 Yellow	DI.851.000 Set of 10 pieces

ID-Color Regular Impression trays with rim, solid

Lower

					
L1 / XL DI.821.010 Purple	L2 / L DI.821.020 Green	L3 / M DI.821.030 Blue	L4 / S DI.821.040 Red	L5 / XS DI.821.050 Yellow	DI.861.060 Set of 6 pieces

Upper

					
U1 / XL DI.831.010 Purple	U2 / L DI.831.020 Green	U3 / M DI.831.030 Blue	U4 / S DI.831.040 Red	U5 / XS DI.831.050 Yellow	DI.861.000 Set of 10 pieces




Made of high quality non-magnetic stainless steel (AISI304).
Micro etched inner surface for better retention



Impression trays regular


Perforated without rim



Lower


No.	Size	A mm	B mm	Code
L1	XXS	55	47	DI.700.010
L2	XS	60	49	DI.700.020
L3	S	72	54	DI.700.030
L4	M	74	57	DI.700.040
L5	L	76	59	DI.700.050
L6	XL	78	61	DI.700.060

Perforated with rim



Lower

No.	Size	A mm	B mm	Code
L1	XXS	55	47	DI.750.010
L2	XS	60	49	DI.750.020
L3	S	72	54	DI.750.030
L4	M	74	57	DI.750.040
L5	L	76	59	DI.750.050
L6	XL	78	61	DI.750.060



Upper

No.	Size	A mm	B mm	Code
U1	XXS	61	47	DI.705.010
U2	XS	64	51	DI.705.020
U3	S	66	53	DI.705.030
U4	M	69	57	DI.705.040
U5	L	71	60	DI.705.050
U6	XL	73	62	DI.705.060



Upper

No.	Size	A mm	B mm	Code
U1	XXS	61	47	DI.755.010
U2	XS	64	51	DI.755.020
U3	S	66	53	DI.755.030
U4	M	69	57	DI.755.040
U5	L	71	60	DI.755.050
U6	XL	73	62	DI.755.060




Sets	Contents	Code
6 pieces	L3,L4,L5,U3,U4,U5	DI.720.060
8 pieces	L3,L4,L5,L6,U3,U4,U5,U6	DI.720.080
12 pieces	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	DI.720.000



Sets	Contents	Code
6 pieces	L3,L4,L5,U3,U4,U5	DI.770.060
8 pieces	L3,L4,L5,L6,U3,U4,U5,U6	DI.770.080
12 pieces	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	DI.770.000

Solid without rim



Lower


No.	Size	A mm	B mm	Code
L1	XXS	55	47	DI.710.010
L2	XS	60	49	DI.710.020
L3	S	72	54	DI.710.030
L4	M	74	57	DI.710.040
L5	L	76	59	DI.710.050
L6	XL	78	61	DI.710.060

Solid with rim




Lower

No.	Size	A mm	B mm	Code
L1	XXS	55	47	DI.760.010
L2	XS	60	49	DI.760.020
L3	S	72	54	DI.760.030
L4	M	74	57	DI.760.040
L5	L	76	59	DI.760.050
L6	XL	78	61	DI.760.060



Upper

No.	Size	A mm	B mm	Code
U1	XXS	61	47	DI.715.010
U2	XS	64	51	DI.715.020
U3	S	66	53	DI.715.030
U4	M	69	57	DI.715.040
U5	L	71	60	DI.715.050
U6	XL	73	62	DI.715.060



Upper

No.	Size	A mm	B mm	Code
U1	XXS	61	47	DI.765.010
U2	XS	64	51	DI.765.020
U3	S	66	53	DI.765.030
U4	M	69	57	DI.765.040
U5	L	71	60	DI.765.050
U6	XL	73	62	DI.765.060



Sets	Contents	Code
6 pieces	L3,L4,L5,U3,U4,U5	DI.725.060
8 pieces	L3,L4,L5,L6,U3,U4,U5,U6	DI.725.080
12 pieces	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	DI.725.000



Sets	Contents	Code
6 pieces	L3,L4,L5,U3,U4,U5	DI.780.060
8 pieces	L3,L4,L5,L6,U3,U4,U5,U6	DI.780.080
12 pieces	L1,L2,L3,L4,L5,L6 U1,U2,U3,U4,U5,U6	DI.780.000




Made of high quality non-magnetic stainless steel (AISI304).
Micro etched inner surface for better retention.




Impression trays regular

Perforated with rim



Lower

No.	Size	A mm	B mm	Code
L1	XL	80	63	DI.800.010
L2	L	78	59	DI.800.020
L3	M	76	57	DI.800.030
L4	S	72	57	DI.800.040
L5	XS	69	50	DI.800.050



Upper

No.	Size	A mm	B mm	Code
U1	XL	81	63	DI.810.010
U2	L	74	60	DI.810.020
U3	M	71	55	DI.810.030
U4	S	67	53	DI.810.040
U5	XS	65	51	DI.810.050


Sets	Contents	Code
6 pieces	L2,L3,L4,U2,U3,U4	DI.850.060
10 pieces	L1,L2,L3,L4,L5 U1,U2,U3,U4,U5	DI.850.000

Solid with rim



Lower

No.	Size	A mm	B mm	Code
L1	XL	80	63	DI.820.010
L2	L	78	59	DI.820.020
L3	M	76	57	DI.820.030
L4	S	72	57	DI.820.040
L5	XS	69	50	DI.820.050




Upper

No.	Size	A mm	B mm	Code
U1	XL	81	63	DI.830.010
U2	L	74	60	DI.830.020
U3	M	71	55	DI.830.030
U4	S	67	53	DI.830.040
U5	XS	65	51	DI.830.050

Sets	Contents	Code
6 pieces	L2,L3,L4,U2,U3,U4	DI.860.060
10 pieces	L1,L2,L3,L4,L5 U1,U2,U3,U4,U5	DI.860.000


Impression trays for children

Perforated with rim



Lower

No.	Size	A mm	B mm	Code
L1	S	54	47	DI.200.010
L2	M	62	51	DI.200.020
L3	L	70	55	DI.200.030



Upper

No.	Size	A mm	B mm	Code
U1	S	55	43	DI.210.010
U2	M	60	45	DI.210.020
U3	L	63	49	DI.210.030


Set	Contents	Code
6 pieces	L1,L2,L3,U1,U2,U3	DI.250.000

Solid with rim



Lower

No.	Size	A mm	B mm	Code
L1	S	54	47	DI.220.010
L2	M	62	51	DI.220.020
L3	L	70	55	DI.220.030



Upper

No.	Size	A mm	B mm	Code
U1	S	55	43	DI.230.010
U2	M	60	45	DI.230.020
U3	L	63	49	DI.230.030

Set	Contents	Code
6 pieces	L1,L2,L3,U1,U2,U3	DI.260.000

Impression trays regular Ehrlicke pattern perforated




Lower

No.	Size	A mm	B mm	Code
UB0	XS	65	46	DI.600.000
UB1	S	67	50	DI.600.010
UB2	M	76	54	DI.600.020
UB3	L	79	56	DI.600.030
UB4	XL	79	61	DI.600.040



DI.640.000
Set of 10 pieces
Contents
UB0,UB1,UB2,UB3,UB4
OB0,OB1,OB2,OB3,OB4



Upper

No.	Size	A mm	B mm	Code
OB0	XS	59	52	DI.610.000
OB1	S	66	52	DI.610.010
OB2	M	71	62	DI.610.020
OB3	L	71	64	DI.610.030
OB4	XL	80	70	DI.610.040




Impression Accessories

Made of high quality non-magnetic stainless steel (AISI304).
Micro etched inner surface for better retention.



Impression trays edentulous

Perforated with rim



Lower


No.	Size	A mm	B mm	Code
L1	XL	84	65	DI.300.010
L2	L	77	64	DI.300.020
L3	M	76	60	DI.300.030
L4	S	68	57	DI.300.040

Upper

No.	Size	A mm	B mm	Code
U1	XL	80	62	DI.310.010
U2	L	77	60	DI.310.020
U3	M	70	60	DI.310.030
U4	S	72	56	DI.310.040

Set	Contents	Code
8 pieces	L1,L2,L3,L4,U1,U2,U3,U4	DI.350.000

Solid with rim



Lower

No.	Size	A mm	B mm	Code
L1	XL	84	65	DI.320.010
L2	L	77	64	DI.320.020
L3	M	76	60	DI.320.030
L4	S	68	57	DI.320.040


Upper

No.	Size	A mm	B mm	Code
U1	XL	80	62	DI.330.010
U2	L	77	60	DI.330.020
U3	M	70	60	DI.330.030
U4	S	72	56	DI.330.040

Set	Contents	Code
8 pieces	L1,L2,L3,L4,U1,U2,U3,U4	DI.360.000

Impression trays anterior depressed

Perforated with rim



Lower


No.	Size	A mm	B mm	Code
L1	S	65	48	DI.500.010
L2	M	72	57	DI.500.020
L3	L	75	57	DI.500.030

Upper

No.	Size	A mm	B mm	Code
U1	S	64	56	DI.510.010
U2	M	67	60	DI.510.020
U3	L	72	63	DI.510.030

Set	Contents	Code
6 pieces	L1,L2,L3,U1,U2,U3	DI.550.000

Solid with rim



Lower

No.	Size	A mm	B mm	Code
L1	S	65	48	DI.520.010
L2	M	72	57	DI.520.020
L3	L	75	57	DI.520.030

Upper

No.	Size	A mm	B mm	Code
U1	S	64	56	DI.530.010
U2	M	67	60	DI.530.020
U3	L	72	63	DI.530.030

Set	Contents	Code
6 pieces	L1,L2,L3,U1,U2,U3	DI.560.000

Impression trays edentulous Ehrlicke pattern perforated



Lower

No.	Size	A mm	B mm	Code
UU1	S	70	48	DI.400.010
UU2	M	70	53	DI.400.020
UU3	L	72	57	DI.400.030

Upper

No.	Size	A mm	B mm	Code
OU1	S	64	46	DI.410.010
OU2	M	68	53	DI.410.020
OU3	L	79	61	DI.410.030



DI.450.000
Set of 6 pieces
Contents
UU1,UU2,UU3, OU1,OU2,OU3



Made of high quality non-magnetic stainless steel (AISI304).
Micro etched inner surface for better retention.



Partial impression trays for crown and bridge work

Perforated with rim



Fig. 30
DI.900.300
Upper left / Lower right



Fig. 31
DI.900.310
Lower left / Upper right



Fig. 32
DI.900.320
Lower jaw



Fig. 33
DI.900.330
Upper jaw



Fig. 99
DI.900.990
Lateral

Solid with rim



Fig. 30
DI.905.300
Upper left / Lower right



Fig. 31
DI.905.310
Lower left / Upper right



Fig. 32
DI.905.320
Lower jaw



Fig. 33
DI.905.330
Upper jaw



Fig. 99
DI.905.990
Lateral

Perforated without rim



Lower left / Upper right

No.	Size	Code
Fig. L2D	S	DI.899.040
Fig. L4D	M	DI.899.050
Fig. L6D	L	DI.899.060



Upper left / Lower right

No.	Size	Code
Fig. L2S	S	DI.899.010
Fig. L4S	M	DI.899.020
Fig. L6S	L	DI.899.030

Perforated with rim



Lower left / Upper right

No.	Size	Code
Fig. L2D	S	DI.901.040
Fig. L4D	M	DI.901.050
Fig. L6D	L	DI.901.060



Upper left / Lower right

No.	Size	Code
Fig. L2S	S	DI.901.010
Fig. L4S	M	DI.901.020
Fig. L6S	L	DI.901.030



Aluminum impression trays regular perforated



Lower

No.	Size	A mm	B mm	Code
L0	XXS	63	44	DI.950.010
L1	XS	69	47	DI.950.020
L2	S	72	48	DI.950.030
L3	M	74	50	DI.950.040
L4	L	75	53	DI.950.050
L5	XL	83	58	DI.950.060



Aluminum partial impression trays for crown & bridge work

Lower left / Upper right



No.	Size	Code
Fig. 1	S	DI.985.010
Fig. 2	M	DI.985.020
Fig. 3	L	DI.985.030



Upper

No.	Size	A mm	B mm	Code
U0	XXS	66	55	DI.960.010
U1	XS	69	56	DI.960.020
U2	S	71	57	DI.960.030
U3	M	73	58	DI.960.040
U4	L	75	62	DI.960.050
U5	XL	78	64	DI.960.060



Upper left / Lower right



No.	Size	Code
Fig. 1	S	DI.980.010
Fig. 2	M	DI.980.020
Fig. 3	L	DI.980.030



DI.970.000
Set of 12 pieces
Contents
L0,L1,L2,L3,L4,L5, U0,U1,U2,U3,U4,U5

D-Color



Silicone buttons for impression trays

Fit any impression tray, easy identification of size



DI.999.000
Assorted pack of 24 buttons
6 colors, 4 of each color

Available in 6 colors



Plastic model bases



Complete upper and lower with hinges

DL.090.010 Standard
DL.090.020 Large

Model and impression shipping cases



DL.093.010
for model

DL.093.020
for impression tray

Rubber model bases

Set of 6 pieces, assorted sizes, upper and lower



DL.091.010