


D16



Photography and X-Ray





Intraoral photographic mirrors


Made of stainless steel






Intraoral photographic mirrors

Double sided front surface
No double reflection
Chromium coated glass
Provides sharp images



- DD.470.000
Set of 7 mirrors and 1 handle:
1 x DD.470.010
1 x DD.470.011
1 x DD.470.020
1 x DD.470.021
1 x DD.470.030
1 x DD.470.031
1 x DD.470.040
1 x DD.472.000



Large mirrors with handle



For intraoral photography
Bendable at any angle
Chrome front surface
No double reflection
Autoclavable

- DD.471.010 Ø 38mm
DD.471.020 Ø 50mm

Handle for photographic mirrors



Designed by
Falcon
DD.472.000
Stainless Steel





Contrastors



DD.475.030
Anterior-Small

DD.475.040
Anterior-Large



DD.475.050
Anterior

DD.475.060
Posterior



DD.475.070
Occlusal-Small



DD.475.080
Occlusal-Medium



DD.475.090
Occlusal-Large

Sterilization cassette for photographic mirrors



DZ.085.090
Suitable for all Intraoral photographic mirrors



Photo cheek retractors



DS.434.010 Small
DS.434.020 Large



DS.435.020
Simplex-Falcon




X-Ray Positioners

Parallelism Technique

- Better Diagnosis, Less Distortion, Avoids Repetitions



Autoclavable



Positioners have special plastic film holders properly aligned to aiming rings. Steel arms are projected to keep them aligned following the Parallelism Technique. Very easy and precise.



X-ray positioners mini kits



DD.1201.00
Anterior



Contents
1 x Indicator arm
1 x Aiming ring
1 x Film holder



DD.1201.02
Aiming ring

DD.1201.01
Film holder



DD.1201.03
Indicator arm



DD.1202.00
Posterior




Contents
1 x Indicator arm
1 x Aiming ring
1 x Film holder



DD.1202.02
Aiming ring

DD.1202.01
Film holder



DD.1202.03
Indicator arm



DD.1203.00
Bite-Wing



Contents
1 x Indicator arm
1 x Aiming ring
1 x Film holder



DD.1203.02
Aiming ring

DD.1203.01
Film holder



DD.1203.03
Indicator arm



X-ray positioners mini kit, Endo



DD.1206.00
Contents
1 x Indicator arm
1 x Aiming ring
1 x Film holder



DD.1206.02
Aiming ring




135° C
Autoclavable

DD.1206.03
Indicator arm



DD.1206.01
Film holder

X-Ray Positioners Set



DD.1200.00
Contents

2 x Indicator arms
2 x Aiming rings
6 x Film holders



2 x Indicator arms
2 x Aiming rings
10 x Film holders



1 x Indicator arm
1 x Aiming ring
4 x Film holders



1 x Indicator arm
1 x Aiming ring
4 x Film holders

X-Ray accessories

Black glasses for X-Ray protection



DZ.9000.04

X-Ray clip



DZ.970.010

X-Ray hangers



Fig. 10
DZ.972.100
10 clips

Fig. 14
DZ.972.140
14 clips

X-Ray protective aprons

A series of products designed to protect against radiation has been designed to combine comfort and uncompromised protection against secondary radiation. Made of lightweight, flexible and multi-lead rubber providing excellent performance.



Type: ST (Standard)



Type: SK (With thyroid collar)



Type: SP (Panoramic)

Lightweight – 0.25 mm Pb on all surface

50 x 80 mm
40 x 60 mm

DD.1360.01
DD.1360.02

DD.1365.01
DD.1365.02

DD.1370.01
DD.1370.02

Heavy – 0.35 mm Pb in front and 0.25 mm on the other parts

50 x 80 mm
40 x 60 mm

DD.1360.03
DD.1360.04

DD.1365.03
DD.1365.04

DD.1370.03
DD.1370.04

Heavy – 0.50 mm Pb in front and 0.25 mm on the other parts

50 x 80 mm
40 x 60 mm

DD.1360.05
DD.1360.06

DD.1365.05
DD.1365.06

DD.1370.05
DD.1370.06

Colors available:



K-01



K-02



K-03



K-04



K-05



K-06




K-07



K-08

Attention: All aprons are made to order. Delivery time 2 weeks.



X-Ray Illuminators

Made in Poland

1 Year Guarantee



PANORAM-01
DQ.8000.01

Screen size: 285 x 155 mm
Dimensions: 358 x 172 x 65 mm
X-Ray film viewer for pantomograms
Desktop or wall mounted



CEFALOM-01
DQ.8010.01

Screen size: 240 x 300 mm
Dimensions: 370 x 330 x 70 mm
X-Ray film viewer for cephalograms
Desktop